

## Recovery Oxygen Mask



Recovery O2 Large 32490B10      Recovery O2 Small #32490B15      Recovery O2 Feline 32490B16  
\*Each mask includes 7 feet of oxygen tubing

The Recovery O2 masks can be used to provide oxygen to animals. The mask has a 22mm x hose barb adapter on the inlet side, allowing you to connect either oxygen tubing (included) or the adapter can be removed to allow the attachment of a standard 22mm breathing circuit. The mask also has dual vents on the sides to allow unrestricted inhalation and exhalation. If oxygen flow is set too low for patient requirement, the vents will allow the patient to inhale ambient air. If oxygen flow is set too high for patient requirement, the vents will allow excess air to escape. The masks and tubing can be cleaned with a dilute chlorhexidine solution or mild soap and water.

Recovery O2 Large #32490B10 (6"H, 5 1/2"D, 2 1/4" hole in diaphragm)

Recovery O2 Small #32490B15 (3 1/2"H, 3 5/8"D, 1 1/4" hole in diaphragm)

Recovery O2 Feline #32490B16 (3 3/4"H, 3 5/8"D, 1 1/8" hole in diaphragm)

Recovery O2 Mask Set V7148 (includes one mask of each size and tubing)